

For Research Use Only

EFEMP1 Polyclonal antibody

Catalog Number: **11150-1-AP**



Basic Information

Catalog Number:

11150-1-AP

Size:

370 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1594

GenBank Accession Number:

BC014410

GeneID (NCBI):

2202

UNIPROT ID:

Q12805

Full Name:

**EGF-containing fibulin-like
extracellular matrix protein 1**

Calculated MW:

55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

