

For Research Use Only

NOTCH2/NOTCH2NL Polyclonal antibody



Catalog Number: 11197-1-AP

1 Publications

Basic Information

Catalog Number:

11197-1-AP

Size:

150 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1679

GenBank Accession Number:

BC019835

GeneID (NCBI):

388677

UNIPROT ID:

PODPK4

Full Name:

Notch homolog 2 (Drosophila) N-terminal like

Calculated MW:

265 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human

Background Information

Notch2 is another homologue gene of Norch in Drosophila. It is a 300kd transmembrane protein, with one 110kd processed fragment, on western blotting. Notch2, like other Norch genes, functions in development of various tissues. This antibody can recognize NOTCH2 and NOTCH2NL.

Notable Publications

Author	Pubmed ID	Journal	Application
Hansen Lin	38462037	Cancer Lett	WB, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

