

For Research Use Only

RBMX Polyclonal antibody

Catalog Number: 11205-1-AP



Basic Information

Catalog Number: 11205-1-AP	GenBank Accession Number: BC007435	Purification Method: Antigen affinity purification
Size: 230 µg/ml	GeneID (NCBI): 27316	
Source: Rabbit	UNIPROT ID: P38159	
Isotype: IgG	Full Name: RNA binding motif protein, X-linked	
Immunogen Catalog Number: AG1659	Calculated MW: 43 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

RBMX, also named as HNRPG, RBMXP1 and Glycoprotein p43, is a rRNA-binding protein which may be involved in pre-mRNA splicing. Investigations into the RBM gene family are uncovering networks of protein interactions which regulate RNA processing, and which might operate downstream of signal transduction pathways.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

