

For Research Use Only

# RBMX Polyclonal antibody

Catalog Number: 11205-1-AP



## Basic Information

<b>Catalog Number:</b> 11205-1-AP	<b>GenBank Accession Number:</b> BC007435	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 230 µg/ml	<b>GeneID (NCBI):</b> 27316	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P38159	
<b>Isotype:</b> IgG	<b>Full Name:</b> RNA binding motif protein, X-linked	
<b>Immunogen Catalog Number:</b> AG1659	<b>Calculated MW:</b> 43 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse, rat

## Background Information

RBMX, also named as HNRPG, RBMXP1 and Glycoprotein p43, is a rRNA-binding protein which may be involved in pre-mRNA splicing. Investigations into the RBM gene family are uncovering networks of protein interactions which regulate RNA processing, and which might operate downstream of signal transduction pathways.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Aliquoting is unnecessary for -20°C storage**

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data