

For Research Use Only

CXCL3/GRO Gamma Polyclonal antibody, PBS Only

Catalog Number: 11221-1-PBS



Basic Information

Catalog Number: 11221-1-PBS	GenBank Accession Number: BC016308	Purification Method: Antigen affinity purification
Concentration: 1 mg/ml	GeneID (NCBI): 2921	
Source: Rabbit	UNIPROT ID: P19876	
Isotype: IgG	Full Name: chemokine (C-X-C motif) ligand 3	
Immunogen Catalog Number: AG1732	Calculated MW: 11 kDa	
	Observed MW: 8-10 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

CXC chemokine ligand 3 (CXCL3) is a member of CXC-type chemokine family that is identified as a major regulator in immune and inflammation responses. CXCL3 is broadly expressed in various human tumor types, and it is also known to play a critical role in mediating tumor development and progression.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Non-treated, LPS and Brefeldin A treated THP-1 cells were subjected to SDS PAGE followed by western blot with 11221-1-AP (CXCL3/GRO Gamma antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11221-1-PBS in a different storage buffer formulation.