

For Research Use Only

PRMT2 Polyclonal antibody

Catalog Number: 11269-1-AP



Basic Information

Catalog Number:

11269-1-AP

Size:

140 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1806

GenBank Accession Number:

BC000727

GeneID (NCBI):

3275

UNIPROT ID:

P55345

Full Name:

protein arginine methyltransferase 2

Calculated MW:

49 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human placenta tissue,

Background Information

PRMT2 (also known as HRMT1L1) is a member of the arginine methyltransferase family which catalyzes arginine methylation and regulates diverse cellular processes, including transcription, translation, DNA repair, and RNA processing. Recently PRMT2 has been identified as a pro-tumorigenic factor in malignant gliomas.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

