

For Research Use Only

NFKBIE Polyclonal antibody

Catalog Number: 11273-1-AP

1 Publications



Basic Information

Catalog Number:

11273-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1800

GenBank Accession Number:

BC011676

GeneID (NCBI):

4794

UNIPROT ID:

O00221

Full Name:

nuclear factor of kappa light

polypeptide gene enhancer in B-cells

inhibitor, epsilon

Calculated MW:

51 kDa

Observed MW:

51 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : human spleen tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yang Zhang	39161335	Oncol Lett	WB, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human spleen tissue were subjected to SDS PAGE followed by western blot with 11273-1-AP (NFKBIE antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.