

TMPRSS4 Polyclonal antibody

Catalog Number: 11283-1-AP

Featured Product

21 Publications

Basic Information

Catalog Number:

11283-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1804

GenBank Accession Number:

BC011703

GeneID (NCBI):

56649

UNIPROT ID:

Q9NRS4

Full Name:

transmembrane protease, serine 4

Calculated MW:

48 kDa

Observed MW:

55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

IHC 1:100-1:400

Applications

Tested Applications:

IHC, WB, ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human, rat, mouse

Cited Species:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: SW 1990 cells, BxPC-3 cells, mouse stomach tissue

IHC: human pancreas cancer tissue, human colon tissue, human colon cancer tissue, human stomach cancer tissue

Background Information

TMPRSS4, also named as CAPH2 and MT-SP2, belongs to the peptidase S1 family. TMPRSS4 is a novel type II transmembrane serine protease that is highly expressed on the cell surface in pancreatic, thyroid and other cancer tissues. It promotes invasion, migration and metastasis of human tumor cells by facilitating an epithelial-mesenchymal transition (EMT) (PMID:20118200).

Notable Publications

Author	Pubmed ID	Journal	Application
Samar Rezq	34575910	Int J Mol Sci	WB
Qin Zhang	33022277	Gastroenterology	WB
Xue-Feng Zhao	34595980	Int J Biol Markers	IHC, WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 4006900926


E: Proteintech-CN@ptglab.comW: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 11283-1-AP (TMPRSS4 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 11283-1-AP (TMPRSS4 Antibody) at dilution of 1:200 (under 40x lens).



SW 1990 cells were subjected to SDS PAGE followed by western blot with 11283-1-AP (TMPRSS4 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.