

For Research Use Only

PLEKHA3 Polyclonal antibody

Catalog Number: 11292-1-AP



Basic Information

Catalog Number:

11292-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1851

GenBank Accession Number:

BC044567

GeneID (NCBI):

65977

UNIPROT ID:

Q9HB20

Full Name:

pleckstrin homology domain
containing, family A
(phosphoinositide binding specific)
member 3

Calculated MW:

34 kDa

Observed MW:

34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : COLO 320 cells, Jurkat cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11292-1-AP (PLEKHA3 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.