

For Research Use Only

# SNAPC1 Polyclonal antibody

Catalog Number: 11297-1-AP



## Basic Information

Catalog Number:

11297-1-AP

Size:

230 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1812

GenBank Accession Number:

BC019038

GeneID (NCBI):

6617

UNIPROT ID:

Q16533

Full Name:

small nuclear RNA activating complex, polypeptide 1, 43kDa

Calculated MW:

43 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human colon cancer using 11297-1-AP (SNAPC1 antibody) at dilution of 1:50 (under 10x lens).