

HLA-DPB1 Polyclonal antibody

Catalog Number: 11298-1-AP

1 Publications

Basic Information

Catalog Number:

11298-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1815

GenBank Accession Number:

BC015000

GeneID (NCBI):

3115

UNIPROT ID:

P04440

Full Name:

major histocompatibility complex,
class II, DP beta 1

Calculated MW:

29 kDa

Observed MW:

29 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

IHC 1:50-1:500

IF-P 1:50-1:500

Applications

Tested Applications:

WB, IHC, IF-P, ELISA

Cited Applications:

IHC

Species Specificity:

human

Cited Species:

rat

Note-IHC: suggested antigen retrieval with
**TE buffer pH 9.0; (*) Alternatively, antigen
retrieval may be performed with citrate
buffer pH 6.0**

Positive Controls:

WB : Daudi cells,

IHC : human tonsillitis tissue, human spleen tissue

IF-P : human tonsillitis tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Karin Harumi Uchima Koecklin	37069178	Sci Rep	IHC

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:


PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage


Selected Validation Data




Daudi cells were subjected to SDS PAGE followed by western blot with 11298-1-AP (HLA-DPB1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.




Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 11298-1-AP (HLA-DPB1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 11298-1-AP (HLA-DPB1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using HLA-DPB1 antibody (11298-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using HLA-DPB1 antibody (11298-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).