

For Research Use Only

# Fas Ligand Polyclonal antibody

Catalog Number: 11322-1-AP

2 Publications



## Basic Information

Catalog Number:

11322-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1816

GenBank Accession Number:

BC017502

GeneID (NCBI):

356

UNIPROT ID:

P48023

Full Name:

Fas ligand (TNF superfamily, member 6)

Calculated MW:

281 aa, 31 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, pig

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Ben-Cheng Lin	31954379	J Biochem Mol Toxicol	WB
Jian Li	38508351	Comp Biochem Physiol C Toxicol Pharmacol	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

