

For Research Use Only

EXT2 Polyclonal antibody, PBS Only

Catalog Number:11348-1-PBS



Basic Information

Catalog Number:

11348-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1748

GenBank Accession Number:

BC010058

GeneID (NCBI):

2132

UNIPROT ID:

Q93063

Full Name:

exostoses (multiple) 2

Calculated MW:

718 aa, 82 kDa

Observed MW:

70-82 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, IF-P, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Exostosin-2 (EXT2) is type II transmembrane glycosyltransferase forming with EXT1 the heterodimeric heparan sulfate polymerase which catalyzes the elongation of the heparan sulfate glycan backbone (PMID: 22660413, 36402845, 36593275).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 11348-1-AP (EXT2 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11348-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human colon cancer using 11348-1-AP (EXT2 antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 11348-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded human colon cancer tissue using EXT2 antibody (11348-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 11348-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded human colon cancer tissue using EXT2 antibody (11348-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 11348-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using EXT2 antibody (11348-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 11348-1-PBS in a different storage buffer formulation.