

For Research Use Only

SRCRB4D Polyclonal antibody

Catalog Number: 11350-1-AP



Basic Information

Catalog Number:

11350-1-AP

Concentration:

350 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1884

GenBank Accession Number:

BC015651

GeneID (NCBI):

136853

UNIPROT ID:

Q8WTU2

Full Name:

scavenger receptor cysteine rich domain containing, group B (4 domains)

Calculated MW:

61 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse kidney tissue, rat kidney tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11350-1-AP (SRCRB4D antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.