

For Research Use Only

# fibrosin Polyclonal antibody, PBS Only

Catalog Number:11356-1-PBS



## Basic Information

Catalog Number:

11356-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1870

GenBank Accession Number:

BC014212

GeneID (NCBI):

64319

UNIPROT ID:

Q9HAH7

Full Name:

fibrosin

Calculated MW:

52 kDa

Observed MW:

48-52 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

Fibrosin is a lymphokine secreted by activated lymphocytes that induces fibroblast proliferation. The FBR1 (or FBS1) gene, encoding fibrosin (or fibrosin-1), was first cloned by Prakash et al. and was found to be able to stimulate several biological activities related to scarring as a fibrogenic lymphokine. Fibrosin is secreted by activated lymphocytes that induce fibroblast proliferation. In vitro application of fibrosin stimulates expression of alpha smooth muscle actin fibers and plays an important role in up regulating the appearance of myofibroblasts during wound healing and possibly in fibrotic diseases. Fibrosin is also implicated in alcohol-induced hepatic fibrosis in vivo.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data