

For Research Use Only

Septin 1 Polyclonal antibody

Catalog Number: 11370-1-AP **1 Publications**



Basic Information

Catalog Number: 11370-1-AP	GenBank Accession Number: BC012161	Purification Method: Antigen affinity purification
Size: 450 µg/ml	GeneID (NCBI): 1731	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q8WYJ6	
Isotype: IgG	Full Name: septin 1	
Immunogen Catalog Number: AG1926	Calculated MW: 367 aa, 42 kDa	
	Observed MW: 40-42 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells,
Cited Applications: WB, IF	
Species Specificity: human, mouse, rat	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Prasad V Phatarpekar	32841357	J Cell Biol	WB, IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 11370-1-AP (Septin 1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.