

For Research Use Only

RPH3A Polyclonal antibody, PBS Only

Catalog Number:11396-1-PBS



Basic Information

Catalog Number:

11396-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1956

GenBank Accession Number:

BC017259

GeneID (NCBI):

22895

UNIPROT ID:

Q9Y2J0

Full Name:

rabphilin 3A homolog (mouse)

Calculated MW:

690 aa, 76 kDa

Observed MW:

76 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

RAB3A is a small G protein that is thought to act in exocytosis of neurotransmitters and hormones in synaptic neurotransmission and cell-cell communication. RPH3A (RPH3A), a RAB3A effector, is a neuronal C2 domain tandem containing protein that is involved in the neuronal transmitter process. Huntington's disease (HD) is a hereditary neurodegenerative disorder characterized by progressive motor disturbances, cognitive defects and dementia. Level of rabphilin 3A is substantially decreased in synapses of most brain regions in R6/1 transgenic mouse model of HD suggesting a role of RPH3A in impaired synaptic transmission.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 11396-1-AP (RPH3A antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 11396-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 11396-1-AP (RPH3A antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 11396-1-PBS in a different storage buffer formulation.



COLO 320 cells were subjected to SDS PAGE followed by western blot with 11396-1-AP (RPH3A antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11396-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed SH-SY5Y cells using RPH3A antibody (11396-1-AP) at dilution of 1:400 and CoraLite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 11396-1-PBS in a different storage buffer formulation.