

For Research Use Only

# STK32C Polyclonal antibody

Catalog Number: 11434-1-AP



## Basic Information

Catalog Number:

11434-1-AP

Concentration:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1960

GenBank Accession Number:

BC015792

GeneID (NCBI):

282974

UNIPROT ID:

Q86UX6

Full Name:

serine/threonine kinase 32C

Calculated MW:

486 aa, 55 kDa

Observed MW:

43-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11434-1-AP (STK32C antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.