

For Research Use Only

STK32C Polyclonal antibody

Catalog Number: 11434-1-AP



Basic Information

Catalog Number: 11434-1-AP	GenBank Accession Number: BC015792	Purification Method: Antigen affinity purification
Size: 133 µg/ml	GeneID (NCBI): 282974	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	UNIPROT ID: Q86UX6	
Isotype: IgG	Full Name: serine/threonine kinase 32C	
Immunogen Catalog Number: AG1960	Calculated MW: 486 aa, 55 kDa	
	Observed MW: 43-55 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse brain tissue, rat brain tissue
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11434-1-AP (STK32C antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.