

For Research Use Only

Flotillin 2 Polyclonal antibody

Catalog Number: 11446-1-AP

2 Publications



Basic Information

Catalog Number:

11446-1-AP

Size:

290 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1996

GenBank Accession Number:

BC017292

GeneID (NCBI):

2319

UNIPROT ID:

Q14254

Full Name:

flotillin 2

Calculated MW:

379 aa, 42 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Flotillins are membrane-associated proteins that contain two homologous isoforms, flotillin 1 (FLOT1) and flotillin 2 (FLOT2). They are highly conserved and widely expressed proteins that commonly used as markers of lipid rafts. Flotillins are involved in various signaling pathways, cell adhesion, membrane trafficking and axonal growth.

Notable Publications

Author	Pubmed ID	Journal	Application
Min Cheng	39265655	Antiviral Res	
Oindrila Paul	38771750	PLoS One	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

