

For Research Use Only

# GPSM3 Polyclonal antibody

Catalog Number: **11470-1-AP**



## Basic Information

Catalog Number:

**11470-1-AP**

Size:

**230 µg/ml**

Source:

**Rabbit**

Isotype:

**IgG**

Immunogen Catalog Number:

**AG2049**

GenBank Accession Number:

**BC018724**

GeneID (NCBI):

**63940**

UNIPROT ID:

**Q9Y4H4**

Full Name:

**G-protein signaling modulator 3  
(AGS3-like, C. elegans)**

Calculated MW:

**160 aa, 18 kDa**

Purification Method:

**Antigen affinity purification**

## Applications

Tested Applications:

**ELISA**

Species Specificity:

**human**

## Background Information

### Storage

Storage:

**Store at -20°C. Stable for one year after shipment.**

Storage Buffer:

**PBS with 0.02% sodium azide and 50% glycerol pH 7.3.**

**Aliquoting is unnecessary for -20°C storage**

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

