

For Research Use Only

EFHC1 Polyclonal antibody

Catalog Number: 11471-1-AP



Basic Information

Catalog Number:

11471-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2048

GenBank Accession Number:

BC020210

GeneID (NCBI):

114327

UNIPROT ID:

Q5JVL4

Full Name:

EF-hand domain (C-terminal)
containing 1

Calculated MW:

640 aa, 74 kDa

Observed MW:

74 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11471-1-AP (EFHC1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.