

For Research Use Only

Cadherin-19 Polyclonal antibody

Catalog Number: 11473-1-AP



Basic Information

Catalog Number:

11473-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2000

GenBank Accession Number:

BC015877

GeneID (NCBI):

28513

UNIPROT ID:

Q9H159

Full Name:

cadherin 19, type 2

Calculated MW:

55 kDa; 87 kDa

Observed MW:

100 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, HeLa cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11473-1-AP (Cadherin-19 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.