

For Research Use Only

ERAL1 Polyclonal antibody

Catalog Number: 11478-1-AP

Featured Product

12 Publications



Basic Information

Catalog Number:

11478-1-AP

Size:

347 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2019

GenBank Accession Number:

BC019094

GeneID (NCBI):

26284

UNIPROT ID:

O75616

Full Name:

Era G-protein-like 1 (E. coli)

Calculated MW:

437 aa, 48 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IF, WB

Species Specificity:

human

Cited Species:

human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Agata Rozanska	28512204	Biochem J	WB
Tian Luo	28553621	Front Cell Infect Microbiol	WB,IF
Di Hu	30877431	Acta Neuropathol	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data