

For Research Use Only

# FARP2 Polyclonal antibody, PBS Only

Catalog Number: 11488-1-PBS



## Basic Information

Catalog Number:

11488-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1990

GenBank Accession Number:

BC021301

GeneID (NCBI):

9855

UNIPROT ID:

O94887

Full Name:

FERM, RhoGEF and pleckstrin domain protein 2

Calculated MW:

120 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data