

For Research Use Only

MED17 Polyclonal antibody

Catalog Number: 11505-1-AP **1 Publications**



Basic Information

Catalog Number: 11505-1-AP	GenBank Accession Number: BC021101	Purification Method: Antigen affinity purification
Size: 400 µg/ml	GeneID (NCBI): 9440	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q9NVC6	
Isotype: IgG	Full Name: mediator complex subunit 17	
Immunogen Catalog Number: AG2057	Calculated MW: 651 aa, 73 kDa	
	Observed MW: 80 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human liver tissue,
Cited Applications: IF	
Species Specificity: human, mouse, rat	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Haining Zhou	26743940	Nat Commun	IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human liver tissue were subjected to SDS PAGE followed by western blot with 11505-1-AP (MED17 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.