

For Research Use Only

MED17 Polyclonal antibody

Catalog Number: 11505-1-AP

1 Publications



Basic Information

Catalog Number:

11505-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2057

GenBank Accession Number:

BC021101

GeneID (NCBI):

9440

UNIPROT ID:

Q9NVC6

Full Name:

mediator complex subunit 17

Calculated MW:

651 aa, 73 kDa

Observed MW:

80 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

IF

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : human liver tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Haining Zhou	26743940	Nat Commun	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human liver tissue were subjected to SDS PAGE followed by western blot with 11505-1-AP (MED17 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.