

For Research Use Only

FNTB Polyclonal antibody

Catalog Number: 11507-1-AP

Featured Product



Basic Information

Catalog Number:

11507-1-AP

Concentration:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2061

GenBank Accession Number:

BC020232

GeneID (NCBI):

2342

UNIPROT ID:

P49356

Full Name:

farnesyltransferase, CAAX box, beta

Calculated MW:

437 aa, 49 kDa

Observed MW:

49 kDa, 44 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

IHC 1:50-1:500

Applications

Tested Applications:

WB, IP, IHC, ELISA

Species Specificity:

human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : mouse heart tissue, A431 cells, mouse skin tissue

IP : HeLa cells,

IHC : human lymphoma tissue,

Background Information

FNTB is the catalytic β subunit of farnesyltransferase which is an alpha/beta heterodimeric enzyme attaching a farnesyl group to a single cysteine in several cellular proteins like RAS and facilitating their translocation from the cytoplasm to the plasma membrane. Farnesyltransferase inhibitors have been developed as ant-tumor agents. FNTB is required for the homeostasis of skin keratinocytes.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 11507-1-AP (FNTB antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



WB result of FNTB antibody (11507-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-FNTB transfected A431 cells.



IP result of anti-FNTB (IP:11507-1-AP, 4ug; Detection:11507-1-AP 1:1000) with HeLa cells lysate 800ug.