

For Research Use Only

PHF15 Polyclonal antibody

Catalog Number: 11513-1-AP **1 Publications**



Basic Information

Catalog Number: 11513-1-AP	GenBank Accession Number: BC021962	Purification Method: Antigen affinity purification
Size: 350 µg/ml	GeneID (NCBI): 23338	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q9NQCC1	
Isotype: IgG	Full Name: PHD finger protein 15	
Immunogen Catalog Number: AG2073	Calculated MW: 791 aa, 88 kDa	
	Observed MW: 110-120 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : MDCK cells,
Cited Applications: WB	
Species Specificity: human, Canine	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Suli Lv	35970393	J Biol Chem	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MDCK cells were subjected to SDS PAGE followed by western blot with 11513-1-AP (PHF15 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.