

For Research Use Only

# ABI1 Polyclonal antibody

Catalog Number: 11516-1-AP



## Basic Information

<b>Catalog Number:</b> 11516-1-AP	<b>GenBank Accession Number:</b> BC024254	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 530 µg/ml	<b>GeneID (NCBI):</b> 10006	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q8IZP0	
<b>Isotype:</b> IgG	<b>Full Name:</b> abl-interactor 1	
<b>Immunogen Catalog Number:</b> AG2077	<b>Calculated MW:</b> 507 aa, 55 kDa	

## Applications

**Tested Applications:**  
IF/ICC, ELISA

**Species Specificity:**  
human, mouse, rat

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of HepG2 cells, using ABL1 antibody 11516-1-AP at 1:100 dilution and PE-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).