

For Research Use Only

# MFSD10 Polyclonal antibody, PBS Only

Catalog Number: 11518-1-PBS



## Basic Information

Catalog Number:

11518-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2080

GenBank Accession Number:

BC014979

GeneID (NCBI):

10227

UNIPROT ID:

Q14728

Full Name:

major facilitator superfamily domain containing 10

Calculated MW:

455 aa, 48 kDa

Observed MW:

48-52 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 11518-1-AP (MFSD10 Antibody) at dilution of 1:600 incubated at 4 degree celsius over night. This data was developed using the same antibody clone with 11518-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 11518-1-AP (MFSD10 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 11518-1-PBS in a different storage buffer formulation.