

For Research Use Only

UCHL5 Polyclonal antibody

Catalog Number: 11527-1-AP

Featured Product

8 Publications



Basic Information

Catalog Number:

11527-1-AP

Concentration:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2095

GenBank Accession Number:

BC025369

GeneID (NCBI):

51377

UNIPROT ID:

Q9Y5K5

Full Name:

ubiquitin carboxyl-terminal
hydrolase L5

Calculated MW:

326 aa, 37 kDa

Observed MW:

38 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

IHC 1:50-1:500

Applications

Tested Applications:

WB, IHC, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

**Note-IHC: suggested antigen retrieval with
TE buffer pH 9.0; (*) Alternatively, antigen
retrieval may be performed with citrate
buffer pH 6.0**

Positive Controls:

WB : mouse colon tissue, mouse brain tissue

IHC : mouse kidney tissue,

Background Information

UCHL5 (Ubiquitin carboxyl-terminal hydrolase isozyme L5) is also named as UCH37 and belongs to the peptidase C12 family. It is functionally linked with multiple protein complexes and signal transduction pathways. UCH37's weight average molecular weight increases as a function of its concentration. By size-exclusion chromatography, the apparent molecular weight of UCH37 was 51, 56, 70, and 109 kDa at 0.15, 0.69, 2.4, and 9.5 mg/mL, respectively. Given UCH37's monomeric molecular weight of 37 kDa and its apparent elongated shape, its tendency toward an apparent molecular weight between 40 and 50 kDa at decreasing concentrations is appropriate for UCH37 monomers (PMID:21953935).

Notable Publications

Author	Pubmed ID	Journal	Application
Po-Ming Chow	36090476	Mol Ther Oncolytics	IHC
Da Liu	32596150	Front Oncol	WB
Zainab Sabry Othman Ahmed	30015387	J Cell Biochem	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse colon tissue were subjected to SDS PAGE followed by western blot with 11527-1-AP (UCHL5 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using 11527-1-AP (UCHL5 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).