

For Research Use Only

CCDC28B Polyclonal antibody

Catalog Number: 11530-1-AP



Basic Information

Catalog Number:

11530-1-AP

Size:

147 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2105

GenBank Accession Number:

BC022848

GeneID (NCBI):

79140

UNIPROT ID:

Q9BUN5

Full Name:

coiled-coil domain containing 28B

Calculated MW:

200 aa, 22 kDa

Observed MW:

30-35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human skeletal muscle tissue, human brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 11530-1-AP (CCDC28B antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.