

For Research Use Only

BRINP1 Polyclonal antibody

Catalog Number: 11544-2-AP **1 Publications**



Basic Information

Catalog Number: 11544-2-AP	GenBank Accession Number: BC021560	Purification Method: Antigen affinity purification
Size: 133 µg/ml	GeneID (NCBI): 1620	
Source: Rabbit	UNIPROT ID: O60477	
Isotype: IgG	Full Name: deleted in bladder cancer 1	
Immunogen Catalog Number: AG2044	Calculated MW: 761 aa, 88 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human

Cited Species:
human

Background Information

BRINP1, also named as DBC1 and FAM5A, inhibits cell proliferation by negative regulation of the G1/S transition. It mediates cell death which is not of the classical apoptotic type and regulates expression of components of the plasminogen pathway.

Notable Publications

Author	Pubmed ID	Journal	Application
Ruibing Chen	28973437	Nucleic Acids Res	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data