

For Research Use Only

GNLY Polyclonal antibody

Catalog Number: 11578-1-AP

2 Publications



Basic Information

Catalog Number:

11578-1-AP

Size:

230 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2163

GenBank Accession Number:

BC023576

GeneID (NCBI):

10578

UNIPROT ID:

P22749

Full Name:

granulysin

Calculated MW:

16 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhichao Wang	34260993	Phytomedicine	WB
Gang Huang	39383633	Phytomedicine	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

