

For Research Use Only

ATPAF2 Polyclonal antibody, PBS Only

Catalog Number: 11598-1-PBS



Basic Information

Catalog Number:

11598-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2164

GenBank Accession Number:

BC032126

GeneID (NCBI):

91647

UNIPROT ID:

Q8N5M1

Full Name:

ATP synthase mitochondrial F1 complex assembly factor 2

Calculated MW:

289 aa, 33 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 11598-1-AP (ATPAF2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11598-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer using 11598-1-AP (ATPAF2 antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 11598-1-PBS in a different storage buffer formulation.