

For Research Use Only

# VPRBP Polyclonal antibody, PBS Only

Catalog Number: 11612-1-PBS

Featured Product



## Basic Information

<b>Catalog Number:</b> 11612-1-PBS	<b>GenBank Accession Number:</b> BC022792	<b>Purification Method:</b> Antigen affinity purification
<b>Concentration:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 9730	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9Y4B6	
<b>Isotype:</b> IgG	<b>Full Name:</b> Vpr (HIV-1) binding protein	
<b>Immunogen Catalog Number:</b> AG2184	<b>Calculated MW:</b> 1506 aa, 169 kDa	
	<b>Observed MW:</b> 169 kDa	

## Applications

**Tested Applications:**  
WB, IHC, IF/ICC, FC (Intra), IP, Indirect ELISA

**Species Specificity:**  
human, mouse

## Background Information

VprBP was first identified as a protein that can interact with HIV-1 viral protein R (PMID: 11223251). It is a component of the CUL4A-RBX1-DDB1-VprBP/DCAF1 E3 ubiquitin-protein ligase complex that could interact with HIV-1 virus Vpr protein and HIV-2 virus Vpx protein (PMID: 18332868; 17314515; 18606781). VprBP is a 1,507-amino acid protein that contains conserved domains, including YXXY repeats, the Lis homology motif, and WD40 repeats. Through binding to Vpr, VprBP allows Vpr to modulate the catalytic activity of the CUL4-DDB1 complex, which in turn leads to the induction of G2 phase arrest in the virus-infected cells (PMID: 17630831). Recently it has been reported that VprBP is able to regulate the p53-induced transcription and apoptotic pathway (PMID: 22184063).

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 11612-1-AP (VPRBP antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11612-1-PBS in a different storage buffer formulation.




WB result of VPRBP antibody (11612-1-AP; 1:8000; incubated at room temperature for 1.5 hours) with sh-Control and sh-VPRBP transfected HEK-293 cells. This data was developed using the same antibody clone with 11612-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 11612-1-AP (VPRBP antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 11612-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 11612-1-AP (VPRBP antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 11612-1-PBS in a different storage buffer formulation.



IP result of anti-VPRBP (IP:11612-1-AP, 3ug; Detection:11612-1-AP 1:1000) with HeLa cells lysate 3800ug. This data was developed using the same antibody clone with 11612-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using VPRBP antibody (11612-1-AP) at dilution of 1:400 and CoraLite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red). This data was developed using the same antibody clone with 11612-1-PBS in a different storage buffer formulation.



$1 \times 10^6$  HeLa cells were intracellularly stained with 0.25 ug VPRBP Polyclonal antibody (11612-1-AP) and CoraLite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2)(red), or 0.25 ug Rabbit IgG control Rabbit PolyAb (30000-0-AP) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 11612-1-PBS in a different storage buffer formulation.