

For Research Use Only

# NRIP Polyclonal antibody

Catalog Number: **11632-1-AP**



## Basic Information

Catalog Number:

**11632-1-AP**

Size:

**130 µg/ml**

Source:

**Rabbit**

Isotype:

**IgG**

Immunogen Catalog Number:

**AG2206**

GenBank Accession Number:

**BC025262**

GeneID (NCBI):

**55827**

UNIPROT ID:

**Q58WW2**

Full Name:

**IQ motif and WD repeats 1**

Calculated MW:

**513 aa, 57 kDa**

Purification Method:

**Antigen affinity purification**

## Applications

Tested Applications:

**WB, ELISA**

Species Specificity:

**human, mouse, rat**

## Background Information

### Storage

Storage:

**Store at -20°C. Stable for one year after shipment.**

Storage Buffer:

**PBS with 0.02% sodium azide and 50% glycerol pH 7.3.**

**Aliquoting is unnecessary for -20°C storage**

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 11632-1-AP (NRIP antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.