

For Research Use Only

NFATC2IP Polyclonal antibody, PBS Only

Catalog Number: 11667-1-PBS



Basic Information

Catalog Number:

11667-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2217

GenBank Accession Number:

BC021551

GeneID (NCBI):

84901

UNIPROT ID:

Q8NCF5

Full Name:

nuclear factor of activated T-cells,
cytoplasmic, calcineurin-dependent 2
interacting protein

Calculated MW:

138 aa, 15 kDa

Observed MW:

50-65 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

NFATC2IP, also named as NIP45, is a 419 amino acid protein, which contains one ubiquitin-like domain. TRAF1 is associated with a fraction of NFATC2IP in the cytoplasm and prevents its translocation to the nucleus. NFATC2IP In T-helper 2 (Th2) cells, regulates the magnitude of NFAT-driven transcription of a specific subset of cytokine genes, including IL3, IL4, IL5 and IL13, but not IL2.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11667-1-AP (NFATC2IP antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11667-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using NFATC2IP antibody (11667-1-AP) at dilution of 1:400 and Multi-Ab CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002), CL594-phalloidin (red). This data was developed using the same antibody clone with 11667-1-PBS in a different storage buffer formulation.