

For Research Use Only

SSX1 Polyclonal antibody

Catalog Number: 11670-1-AP



Basic Information

Catalog Number:

11670-1-AP

Size:

850 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2254

GenBank Accession Number:

BC001003

GeneID (NCBI):

6756

UNIPROT ID:

Q16384

Full Name:

synovial sarcoma, X breakpoint 1

Calculated MW:

22 kDa

Observed MW:

21-25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : PC-3 cells, HeLa cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 11670-1-AP (SSX1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.