

For Research Use Only

CARHSP1 Polyclonal antibody, PBS Only

Catalog Number: 11672-1-PBS

Featured Product



Basic Information

Catalog Number:

11672-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2265

GenBank Accession Number:

BC003366

GeneID (NCBI):

23589

UNIPROT ID:

Q9Y2V2

Full Name:

calcium regulated heat stable protein
1, 24kDa

Calculated MW:

147 aa, 16 kDa

Observed MW:

~22 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

CARHSP1 (calcium-regulated heat-stable protein 1, a 24 kDa protein also known as CRHSP-24) is a ubiquitously expressed phosphoprotein that is comprised of 147 amino acids with nearly 14% proline (PMID: 21990353, 30257675). CARHSP1 was discovered as a biomarker for diabetic retinopathy and as a regulator of TNF- α mRNA stability in macrophages (PMID: 21078874). CARHSP1 may be involved in oxidative stress via traffic between stress granules and processing bodies (PMID: 21177848).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11672-1-AP (CARHSP1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11672-1-PBS in a different storage buffer formulation.



WB result of CARHSP1 antibody (11672-1-AP; 1:8000; incubated at room temperature for 1.5 hours) with sh-Control and sh-CARHSP1 transfected HEK-293 cells. This data was developed using the same antibody clone with 11672-1-PBS in a different storage buffer formulation.