

For Research Use Only

SRP54 Polyclonal antibody

Catalog Number: 11729-1-AP **5 Publications**



Basic Information

Catalog Number: 11729-1-AP	GenBank Accession Number: BC003389	Purification Method: Antigen affinity purification
Size: 387 µg/ml	GeneID (NCBI): 6729	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: P61011	IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate
Isotype: IgG	Full Name: signal recognition particle 54kDa	IHC 1:50-1:500
Immunogen Catalog Number: AG2327	Calculated MW: 54 kDa	
	Observed MW: 54 kDa	

Applications

Tested Applications: IHC, IP, WB, ELISA	Positive Controls:
Cited Applications: WB	WB : human placenta tissue, NIH/3T3 cells
Species Specificity: human, mouse	IP : human placenta tissue,
Cited Species: human	IHC : human breast cancer tissue, human placenta tissue
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

The signal recognition particle (SRP) is a ribonucleoprotein complex that mediates the targeting of proteins to the endoplasmic reticulum (ER). The complex consists of a 7S (or 7SL) RNA and 6 different proteins, and signal recognition particle 54 (SRP54) is one of them. SRP54 binds to the signal sequence of presecretory protein as they emerge from the translating ribosomes, and then transfers them to translocating chain-associating membrane protein (TRAM).

Notable Publications

Author	Pubmed ID	Journal	Application
Dianrong Li	31420216	Mol Cell	
Masashi Ueda	24388981	Biochem Biophys Res Commun	WB
Ren Yong-Feng YF	23468895	PLoS One	WB


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 11729-1-AP (SRP54 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



IP result of anti-SRP54 (IP:11729-1-AP, 4ug; Detection:11729-1-AP 1:500) with human placenta tissue lysate 4000ug.



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 11729-1-AP (SRP54 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).