

For Research Use Only

RPP40 Polyclonal antibody

Catalog Number: 11736-1-AP



Basic Information

Catalog Number:

11736-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2341

GenBank Accession Number:

BC017871

GeneID (NCBI):

10799

UNIPROT ID:

O75818

Full Name:

ribonuclease P/MRP 40kDa subunit

Calculated MW:

35 kDa, 41 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:5000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, COLO 320 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 11736-1-AP (RPP40 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.