

For Research Use Only

DCLRE1C Polyclonal antibody

Catalog Number: 11747-1-AP



Basic Information

Catalog Number:

11747-1-AP

Size:

167 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2343

GenBank Accession Number:

BC022254

GeneID (NCBI):

64421

UNIPROT ID:

Q96SD1

Full Name:

DNA cross-link repair 1C (PSO2 homolog, *S. cerevisiae*)

Calculated MW:

692 aa, 78 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, human kidney tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 11747-1-AP (DCLRE1C antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.