

For Research Use Only

HFE2 Polyclonal antibody

Catalog Number: 11758-1-AP



Basic Information

Catalog Number: 11758-1-AP	GenBank Accession Number: BC017926	Purification Method: Antigen affinity purification
Size: 550 µg/ml	GeneID (NCBI): 148738	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	UNIPROT ID: Q6ZVN8	
Isotype: IgG	Full Name: hemochromatosis type 2 (juvenile)	
Immunogen Catalog Number: AG2354	Calculated MW: 426 aa, 45 kDa	
	Observed MW: 23 kDa, 38 kDa, 50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue,
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 11758-1-AP (HFE2 Antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours.