

For Research Use Only

HFE2 Polyclonal antibody

Catalog Number: 11758-1-AP

1 Publications



Basic Information

Catalog Number:

11758-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2354

GenBank Accession Number:

BC017926

GeneID (NCBI):

148738

UNIPROT ID:

Q6ZVN8

Full Name:

hemochromatosis type 2 (juvenile)

Calculated MW:

426 aa, 45 kDa

Observed MW:

23 kDa, 38 kDa, 50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : mouse liver tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yong Zuo	38218852	Cell Death Dis	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 11758-1-AP (HFE2 Antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours.