

For Research Use Only

# FCN1 Polyclonal antibody

Catalog Number: 11775-1-AP

3 Publications



## Basic Information

Catalog Number:

11775-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2444

GenBank Accession Number:

BC020635

GeneID (NCBI):

2219

UNIPROT ID:

O00602

Full Name:

ficolin (collagen/fibrinogen domain containing) 1

Calculated MW:

326 aa, 35 kDa

Observed MW:

33-37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : human plasma, K-562 cells, human liver tissue, mouse liver tissue, rat liver tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Muhammad Khan	36389725	Front Immunol	IHC
Kohei Kusakari	36052069	Front Immunol	WB
Fangfang Zhong	39139822	Front Genet	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11775-1-AP (FCN1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.