

For Research Use Only

AFTPH Polyclonal antibody

Catalog Number: **11785-1-AP**



Basic Information

Catalog Number:

11785-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2402

GenBank Accession Number:

BC022247

GeneID (NCBI):

54812

UNIPROT ID:

Q6ULP2

Full Name:

aftiphilin

Calculated MW:

937 aa, 102 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Aftiphilin is a component of the clathrin machinery in neurons, binding to clathrin adaptor protein 1 and 2. Western assay showed that the protein emigrate as a 126 kDa protein in SDS PAGE, with larger molecular weight than the predicted size (99 and 102 kDa for the two isoforms).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

