

For Research Use Only

DIO1 Polyclonal antibody, PBS Only

Catalog Number:11790-1-PBS



Basic Information

| | | |
|--|--|--|
| Catalog Number: 11790-1-PBS | GenBank Accession Number: BC017955 | Purification Method: Antigen affinity purification |
| Size: 1 mg/ml | GeneID (NCBI): 1733 | |
| Source: Rabbit | UNIPROT ID: P49895 | |
| Isotype: IgG | Full Name: deiodinase, iodothyronine, type I | |
| Immunogen Catalog Number: AG2410 | Calculated MW: 249 aa, 29 kDa | |
| | Observed MW: 29 kDa, 13 kDa | |

Applications

Tested Applications:
WB, IHC, IF/ICC, FC (Intra), Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

DIO1(Type I iodothyronine deiodinase) is also named as ITDI1, TXDI1,5DI,DIOI and belongs to the iodothyronine deiodinase family.It is a membrane selenoprotein that catalyzes the deiodination of L-thyroxine to the biologically active thyroid hormone 3,3',5-triiodothyronine.DIO1 is located in liver, kidney, and thyroid,which has both ORD and IRD activities(PMID: 9389495). It has some isoforms produced by alternative splicing with the MW of 29 kDa, 21 kDa, 13 kDa, 23 kDa, 19 kDa, 4 kDa, 13 kDa, 9 kDa and 10 kDa.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




mouse lung tissue were subjected to SDS PAGE followed by western blot with 11790-1-AP (DIO1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 11790-1-PBS in a different storage buffer formulation.




Immunofluorescent analysis of (-20°C Ethanol) fixed A549 cells using DIO1 antibody (11790-1-AP) at dilution of 1:200 and CoraLite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 11790-1-PBS in a different storage buffer formulation.



1×10^6 HeLa cells were intracellularly stained with 0.25 ug DIO1 Polyclonal antibody (11790-1-AP) and CoraLite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2)(red), or 0.25 ug Rabbit IgG control Rabbit PolyAb (30000-0-AP) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 11790-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 11790-1-AP (DIO1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 11790-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 11790-1-AP (DIO1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 11790-1-PBS in a different storage buffer formulation.