

For Research Use Only

RGS18 Polyclonal antibody

Catalog Number: 11866-1-AP **1 Publications**



Basic Information

Catalog Number: 11866-1-AP	GenBank Accession Number: BC020632	Purification Method: Antigen affinity purification
Size: 240 µg/ml	GeneID (NCBI): 64407	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9NS28	
Isotype: IgG	Full Name: regulator of G-protein signaling 18	
Immunogen Catalog Number: AG2445	Calculated MW: 235 aa, 28 kDa	
	Observed MW: 28 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human peripheral blood platelets,
Cited Applications: WB	
Species Specificity: human	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yong Meng	38214344	Int J Mol Med	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human peripheral blood leukocyte sample were subjected to SDS PAGE followed by western blot with 11866-1-AP (RGS18 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.